


VERSAFLO3[®] MATTRESS SYSTEM

Feel Rating	Variable
Mattress Height	200mm
Risk Category	Very High Risk
User Weight	20kg - 140kg
Warranty	2 Year Warranty



BED CODES & SIZES

LONG SINGLE

Code: ICVF3LS
Mattress Size: 2030 x 900

KING SINGLE

Code: ICVF3KS
Mattress Size: 2030 x 900

LONG DOUBLE

Code: ICVF3D
Mattress Size: 2030 x 1350

QUEEN

Code: ICVF3Q
Mattress Size: 2030 x 1520

KEY FEATURES

- 200mm Mattress Height
- 140kg Weight Rating
- Anti slip base layer
- Easy clean ZeroTec cover with heavy duty zipper
- CPR emergency release
- Flip free design
- Designed for use on adjustable beds with Contour Latch technology
- 2 Year Warranty

SPECIFICATIONS

- Hybrid alternating air system
- ActiveX layer for comfort and pressure care
- XP Support anti bottom out layer
- Variable pressure and comfort settings
- Variable interval cell pressure settings
- Anti rip base fabric with contour latches
- 20 - 140kg user weight
- Medical Cover MVTR (g/24h·m²): 1635.0 Average
- Medical Cover BS7175: Crib 5
- ARTG Identifier: 478205
- Waterlow ≤ 20

Circulation Intervals	5min, 10min, 15min, 20min, static
Air Capacity	1.5~8L

PRESSURE CARE

Risk Category Very High

This rating is based on our average peak pressure of <16.5mmHg.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

CARE & CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 65°C for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 60°C. The foam can be autoclaved at 134°C.

VersaFlo3[®] comes with repair kit, can be repaired up to five times before requiring replacement.